



Easy® Sentinel™ is the answer to your nurse monitoring needs in the LTM/EMU environment. This solution allows for a staff member to easily monitor up to six studies remotely that are being acquired on the network with Cadwell® Easy III equipment.



Remotely monitor live streaming video and patient event notifications while Easy III records and trends EEG data

Sentinel Easily Integrates into Existing Easy III Environments

Easy Sentinel nurse monitoring stations are fully compatible with the Easy III product line. This station can monitor patients in the ICU on Envisio EEG systems, routine EEG tests occurring in the neurology department, and patients in the EMU. The systems simply need to be running Easy III software with video and be connected to the network. Cadwell offers a variety of options when it comes to Easy III, from ambulatory systems all the way up to a full EMU solution.



Configurable Audible and Visual Notifications Settings

Audible and visual notifications are available and can be associated with a patient button press. Additionally, there is an option to notify the station when spike and seizure events are detected.

The .wav files used for the audible notification can be customized, simply browse to the .wav file and apply the corresponding notification. For example, each patient room could have a custom audio file that clearly states the room number.



Simplified Workflow for Staff

Easy III software is easy to use. Simply double click the icon on the computer desktop, then select the active recordings occurring on the network that need to be monitored. It's easy to add recordings once monitoring. Simply click the add recording button and select the additional patient(s) to monitor.

For more information on how Easy Sentinel will work for you, please give us a call at 800.245.3001.

909 N. Kellogg St. · Kennewick, WA 99336
(800) 245-3001 · (509) 735-6481 ph · (509) 783-6503 fx
www.cadwell.com · info@cadwell.com

