


Technomed Needles Special Pricing



Technomed Disposable Subdermal Needle Electrodes **Now \$35.00***

Standard single stainless steel needle. Multi-color (five) PVC lead wires that terminate into a safety DIN 42 802 connector. Sold in boxes containing five of each color (25 total electrodes). Pre-sterilized using ETO gas. **Single patient use only.**

Length		Diameter		Wire Length		Qty	Order Number
mm	in	mm	Gauge	meters	feet		
13	0.5	0.40	27	2.5	8	1 box (25)	202303-000

*This special applies to the 2.5 meter length only.

Call 800-245-3001 to place your order today!
Visit our website www.cadwell.com to check out our specials each month!!



Advancing medical technology to help you, help others.

06/15/09