

Apollo™ EEG



Apollo[™] EEG

REMOTE MONITORING SOLUTION

Apollo EEG allows physicians and technicians to remotely view and monitor the quality of ambulatory EEG data with HD video.

Compact and only 1 lb. 3 oz., Apollo is designed for patient comfort. The Apollo recorder, the 32-channel Apollo amplifier, the Q-Video[®] Mobile 3 camera with automatic infrared lighting and the NUC PC complete the system.

Apollo provides up to 96 hours of uninterrupted recording on a single charge and allows patients to move freely during continuous data acquisition.



Cadwell Provides

- Arc Apollo Ambulatory System
 - Apollo Recorder and USB cable
 - Apollo Amplifier and Apollo cable
 - Q-Video Mobile 3, USB cable and power cable
- NUC PC and power cable
- Carrying case

Customer Provides

- Mobile Hotspot and power cable
- Remote Access Software (e.g. ScreenConnect or LogMeIn)