



Disposable MR Conditional / CT-Compatible EEG Cup Electrodes

Cadwell® MR Conditional / CT-compatible disposable EEG cup electrodes enable faster resumption of studies after imaging, saving time during long-term epilepsy monitoring and critical care continuous EEG.

- FDA cleared as MR Conditional for 1.5T and 3.0T.
- Radiolucent materials minimize CT artifacts.
- Electrodes remain in place during imaging, saving up to an hour per imaging session.
- Available as both time-saving arrays and singles.



Color-Coded

After imaging, clinical staff can easily reconnect using color-coded labels.

Non-MR Label

Detach the clearly marked non-MR portion of the lead before MRI or CT scans.

Shorter Monitoring Gaps

Reduces non-reimbursable time and avoids long gaps where a seizure could be missed.


Product and Ordering Details

Single Electrodes — MR Conditional / CT-compatible

Available in three electrode configurations and four lead lengths.

**MR Conditional / CT-compatible portion of lead wire is 28 cm (11")*

Part #	Configuration	Lead Length
302932-000-039	28 electrodes	1.0 m (39")
302932-000-059	28 electrodes	1.5 m (59")
302932-000-078	28 electrodes	2.0 m (78")
302932-000-098	28 electrodes	2.5 m (98")
302933-000-039	24 electrodes	1.0 m (39")
302933-000-059	24 electrodes	1.5 m (59")
302933-000-078	24 electrodes	2.0 m (78")
302933-000-098	24 electrodes	2.5 m (98")
302934-000-039	12 electrodes	1.0 m (39")
302934-000-059	12 electrodes	1.5 m (59")
302934-000-078	12 electrodes	2.0 m (78")
302934-000-098	12 electrodes	2.5 m (98")



Arrays — MR Conditional / CT-compatible

Available in two electrode configurations and two lead lengths.

**MR Conditional / CT-compatible portion of array wire is 18 cm (7")*

Part #	Configuration	Lead Length
302930-000-039	3 arrays, total of 27 electrodes	1.0 m (39")
302930-000-059	3 arrays, total of 27 electrodes	1.5 m (59")
302931-000-039	3 arrays, total of 25 electrodes	1.0 m (39")
302931-000-059	3 arrays, total of 25 electrodes	1.5 m (59")



ORDER ONLINE

estore.cadwell.com

CONTACT SUPPLIES & ACCESSORIES

1.800.245.3001 • +1.509.735.6481 • orders@cadwell.com



© 2026 Cadwell Industries, Inc. All rights reserved. Information and products displayed in this document demonstrate sample configurations containing optional components, which may be changed without notice and do not define what is delivered with an order. Actual product configuration and content are determined and confirmed independently at the time of purchase. Product and services availability may vary by country or market. Please contact Cadwell for additional information. This document contains trademarks that belong to Cadwell Industries, Inc. and other companies, respectively. **PN #316011-000 Rev. 00**